



Vitrex® Micro Pipettes (end-to-end)

VITREX end-to-end pipettes are for users who require easy and very accurate capillary tube volume measurement. The pipettes are manufactured to ISO 7550.

Packaging

Glass vials containing 100 pcs. (50 µl - 50 pcs.)
The vials have a dispenser top, making it easy to shake out one tube at a time, thereby ensuring hygienic handling.

Capillary Tube Dispenser
Alcohol Swabs

See page 14
See page 16

The dispenser tops are colour-coded to indicate the tube volume in conformity with ISO 1769. Ten vials are packaged in a sturdy carton in order to protect the tubes from damage during handling or transport.

VITREX tolerances

	Acc	CV
0.5 - 10 µl	± 0.50%	0.75%
20 - 100 µl	± 0.50%	0.60%



Technical Specs.

Vol. (µl)	Colour code	pcs	L (mm) +/- 1mm	Plain	Natrium Hep.	Ammon. Hep.	EDTA Hep.
0.5	○	10x100	15	177613	*	*	*
1	●	10x100	29	165013	167513	169013	*
2	●	10x100	29	165113	167613	169113	*
3	●	10x100	29	165213	167713	169213	*
4	●	10x100	29	165313	167813	169313	*
5	○	10x100	29	165413	167913	169413	170113
10	●	10x100	29	174313	177113	179213	174113
20	●	10x100	29	174513	177213	179313	174413
25	○	10x100	29	174613	177313	179413	178213
40	●	10x100	29	174813	177413	179513	*
44.7	●	10x50	29	175013	177713	178613	*
50	●	10x50	29	174913	177513	179613	171113

* available upon request
For lithium hep. please contact us

Accessories

	Article no.
Pipette holder	140313

Shelf life 4 years

